



System Requirements

STANDALONE & CLIENT SYSTEM REQUIREMENTS		
	Minimum Requirements	Recommended Requirements
Operating System	Windows 7 64-bit / Mac OS X 10.10	Windows 10 64-bit / Mac OS X 10.12
Network (LAN)	100BaseT	1000BaseT
Network (WAN) for Client	10 Mbps down / 2.5 Mbps Up	1000 Mbps down / 100 Mbps Up
Display	1280 x 800 resolution	1920 X 1080 resolution or higher
Processor	Intel Core i5	Intel i7 / Xeon at 2.2 GHZ or higher
Video Card	Nvidia or AMD dedicated GPU	Nvidia dedicated GPU
System RAM	16GB	16GB
Video RAM	2GB	4GB
SYSTEMS WITHOUT DEDICATED GRAPHICS		
Processor / Graphics	Intel i7 with HD 5000 Integrated Graphics	

SERVER REQUIREMENTS		
	Minimum Requirements	Recommended Requirements
Operating System	Windows Server 2008 R2 / Mac OS X 10.10	Windows Server 2012 R2 / Mac OS X 10.12
Network (LAN)	100BaseT	1000BaseT
Network (WAN) for Client	5 Mbps down / 25 Mbps Up	100 Mbps down / 1000 Mbps Up
System RAM	8GB	16GB
Storage Space	1 TB hard drive	2 TB Raid 1+0,5,6, or SAN

About Circle

Circle Cardiovascular Imaging develops and markets cardiovascular post-processing software that allows for the evaluation and analysis of MR and CT images. Available for clinical and research use, the stand-alone software provides full DICOM and PACS connectivity.

Visit circlevi.com for a free trial