

System Requirements

| cvi42 STANDALONE & CLIENT SYSTEM REQUIREMENTS | | | | |
|---|---|-----------|--------------------------------------|-----------|
| | Minimum Requirements | | Recommended Requirements | |
| Operating System | Windows 7 64-bit / Mac OS X 10.10 | | Windows 10 64-bit / Mac OS X 10.12 | |
| Network (LAN) | 100BaseT | | 1000BaseT | |
| Network (WAN) for Client | 10 Mbps down / 2.5 Mbps Up | | 1000 Mbps down / 100 Mbps Up | |
| Display | 1280 x 800 resolution | | 1920 X 1080 resolution or higher | |
| Processor | Intel Core i5 | | Intel i7 / Xeon at 2.2 GHZ or higher | |
| Video Card | Nvidia or AMD dedicated GPU | | Nvidia dedicated GPU | |
| | MR | CT | MR | CT |
| System RAM | 4GB | 8GB | 8GB | 16GB |
| Video RAM | 1GB | 1GB | 2GB | 4GB |
| SYSTEMS WITHOUT DEDICATED GRAPHICS | | | | |
| Processor / Graphics | Intel i7 with HD 3000 Integrated Graphics | | | |

| SERVER REQUIREMENTS | | | | |
|--------------------------|---|-----------|---|-----------|
| | Minimum Requirements | | Recommended Requirements | |
| Operating System | Windows Server 2008 R2 / Mac OS X 10.10 | | Windows Server 2012 R2 / Mac OS X 10.12 | |
| Network (LAN) | 100BaseT | | 1000BaseT | |
| Network (WAN) for Client | 5 Mbps down / 25 Mbps Up | | 100 Mbps down / 1000 Mbps Up | |
| Processor | Intel Core i5 | | Intel i7 / Xeon at 2.2 GHZ or higher | |
| | MR | CT | MR | CT |
| System RAM | 4GB | 8GB | 8GB | 16GB |
| Storage Space | 1 TB hard drive | | 2 TB Raid 1+0,5,6, or SAN | |

About Circle

Circle Cardiovascular Imaging develops and markets cardiovascular post-processing software that allows for the evaluation and analysis of MR and CT images. Available for clinical and research use, the stand-alone software provides full DICOM and PACS connectivity.

Visit circlevi.com for a free trial