


CVI42[®]

Minimum Specifications

Standalone / Client



CPU

Intel Core i5
7th gen +



RAM

8 GB +



VIDEO

Intel Integrated
HD Graphics 630+



DISK

60 MB/s
SSD/NVMe



MAC OS

10.13 High Sierra
10.14 Mojave



WINDOWS


Windows 8.1
Windows 10



NETWORK SPEED

1
Gb/s

Server



CPU

Intel Core i7 7th gen+
or Xeon E5-26xx+



RAM

16 GB



DISK

60 MB/s
RAID or SSD/NVMe



WINDOWS

Server 2012 or 2012 R2
Server 2016
Server 2019



NETWORK SPEED

10
Gb/s

CT & 4D Flow

Higher Requirements



RAM

16 GB+



VIDEO

Discrete nVidia GeForce900 series,
or GeForce 1000 or 2000 series (2GB+ Ram)
or AMD Radeon Rx300 series+ (2GB+ RAM)



USERS

up to 3
Concurrent

Please contact your
Circle advisor for
client-server
configurations
above 3
concurrent users

Questions?

✉ support@circlecvi.com

☎ +1 (403) 338 1870